

High Performance Affinity Chromatography Column

TSKgel® FcR-IIIA-NPR



The First-in-Class, rhFcyreceptor IIIA immobilized affinity chromatography column for separation/analysis of antibody based on the difference of ADCC activity.

Features

- Innovative high performance affinity chromatography column for analysis of antibody
- Recombinant human Fcyreceptor IIIA is immobilized on the surface of non-porous resin as ligand
- Discriminate conformational change of Fc region of antibody attributed by N-linked glycan, which enables to separate antibodies based on the difference of ADCC activity
- Applicable to intact antibodies without any time-consuming or burdensome pretreatment, purification
- Rapid separation of antibody in 30 min is available

Product

P/N	Description	Particle size	dimension
0023513	TSKgel FcR-IIIA-NPR	5 µm	4.6 mm I.D. × 7.5 cm

* Guard column is not available. * Recommend to install line filter between pump and injection valve, Line filter kit PEEK, P/N 0018014.


Target molecules and applications

- Antibody ○ Fc fusion protein
- Screening of cell line ● Optimization of culturing condition
- Monitoring of upstream process ● Quality control of mAb

Basic properties


Description	TSKgel FcR-IIIA-NPR
Base materials	Hydrophilic non-porous polymer
Particle size	5 µm
Ligand	Recombinant human Fcyreceptor IIIA
Column size	4.6 mm I.D. × 7.5 cm
Column hardware	PEEK

Concept



MAb therapeutics

N-glycan structures in Fc region would affect drug efficacy, and might cause lot-to-lot variation




Chromatogram

TSKgel FcR-IIIA-NPR discriminates glycan structure and separates IgG based on ADCC activity

Applications


• Separation of mAbs



Conditions 1
 Column: TSKgel FcR-III-A-NPR
 (4.6 mm I.D. × 7.5 cm)
 Eluent A: 50 mmol/L sodium citrate buffer (pH 6.5)
 B: 50 mmol/L sodium citrate buffer (pH 4.5)
 Gradient: B 0 – 100 %
 (2 – 20 min, linear)
 Flow rate: 1.0 mL/min
 Detection: UV (280 nm)
 Temp.: 25 °C
 Samples: monoclonal antibody

Distinctive chromatograms were obtained for all monoclonal antibodies applied on TSKgel FcR-III-A-NPR column. (Typical monoclonal antibody therapeutics produced by CHO cell shows 3 peaks)


• Lot-to-lot variability of packing materials



Conditions 2
 Conditions are the same as conditions 1

Almost no difference in chromatogram was observed in gel lots


• Column durability (Separation of antibody in cell culture)



Conditions 3
 Eluent A: 50 mmol/L sodium citrate buffer + 150 mmol/L NaCl (pH 6.5)
 B: 50 mmol/L sodium citrate buffer + 150 mmol/L NaCl (pH 4.5)
 Gradient: B 0 – 100 % (7 – 25 min, linear)
 Temp.: 20 °C
 Samples: CHO cell culture supernatant (1 g/L, mAb), purified mAb
 *Sample was filtered and applied
 *Other conditions were the same as conditions 1

Applicable to mAb in CHO cell culture with good reproducibility


• Application to quality control (Monitoring of lot-to-lot variability of mAb therapeutics)



Conditions 4
 Samples: Commercial mAb therapeutics (2 lots)
 * Other conditions were the same as conditions 1

Applicable to lot-to-lot comparison with ease and rapidly

• Application to process analysis (Analysis of mAb cultured in CHO cell)



Conditions 5
 Sample: mAb, 5 µg
 * CHO cell culture was sampled periodically after starting culture, filtered, purified by protein A resin and applied on TSKgel FcR-III-A-NPR
 * Other conditions were the same as conditions 1
 * CHO cell culture was kindly provided by Manufacturing Technology Association of Biologics

Applicable to monitor the time course of mAb during culture



TOSOH

* TSKgel, TOYOPEARL are registered trademarks of Tosoh Corporation in Japan, United States, European Union, etc.
 * Specifications and design may be changed without prior notice.

This research is partially supported by the developing key technologies for discovering and manufacturing pharmaceuticals used for next-generation treatments and diagnoses from Japan Agency for Medical Research and Development (AMED).
 Grant Number JP17ae0101003.

TOSOH CORPORATION

Bioscience Division

3-8-2 Shiba, Minato-ku, Tokyo 105-8623, Japan

Phone : +81 3 5427 5180

Fax : +81 3 5427 5220

This brochure may not be reprinted or copied in whole or in part without written consent from Tosoh Corporation.
 Copyright 2018 by Tosoh Corporation