

Vaccine purity for safety and efficacy

Validated and ready-to-use Host Cell Protein ELISA Kits



Our most popular HCP ELISA Kits

• HeLa (new)	• MDCK (new)
• HEK 293	• PER.C6®
• A549	• CAP®
• Vero	• BHK
• SF9	• MRC5

Trusting your impurity analysis is essential for vaccine purity characterization. Cygnus Technologies HCP ELISA kits deliver the exceptional reliability you need. With over 35 impurity analysis ELISAs, bioprocessing impurity analysis is what we do. If you don't see the assay you need, we can create one for you.

- **Reliable and reproducible:** <10% CV
- **Compliant:** validated per regulatory guidelines
- **Highly sensitive:** as low as 0.2 ng/mL LLOQ
- **Comprehensive:** detect all HCP analytes in a single assay

Improve your vaccine safety and efficacy. Trust the leader in bioprocess impurity analysis. Find out more at cygnustechnologies.com/hcp

PER.C6 is a registered trademark of Crucell Holland B.V.
CAP is a registered trademark of CEVEC Pharmaceuticals GmbH

© 2017 Cygnus Technologies. All rights reserved.

For research use only. Not intended for animal or human therapeutic or diagnostic use.