



MEDIUM KULTUR SEL PUNCA & SEL DEWASA



Miltenyi Biotec

Empowering cell research. Enabling cell therapy.

Xeno & Serum-Free System Media For Human MSCs, HSC, NK Cells & T Cells.

- High-quality cell culture media for consistency and accuracy
- Manufactured according to strict quality standards
- Validated stem cell products and growth factors
- Xeno & Serum - Free





Produk medium Miltenyi Biotec (Germany) dapat digunakan untuk melakukan kultur sel punca berupa MSC, HSC, serta sel-sel NK dan sel-sel T. Medium tersebut dilengkapi dengan supplement yang mengandung *Growth Factor* serta tambahan lainnya.



Medium	StemMACS MSC Expansion Media Kit XF	StemMACS HSC Expansion Cocktail	NK MACS Medium	TexMACS Medium
Reactivity	Human Mesenchymal Stem Cells (MSCs)	CD34 ⁺ hematopoietic stem and progenitor cells	Human NK Cells	Human and Mouse T cells and regulatory T cells.
Source	Bone marrow, Lipoaspirate and etc	Cord blood, Peripheral Blood, Bone Marrow, etc	Human Blood	Human Blood
Size	500 mL	500 mL & 100 mL	500 mL	500 mL
Supplement	+ [7 mL]	-	+ [5 mL]	
Xeno & Serum Free	+	+	without animal-derived components without antibiotics	without animal-derived components without antibiotics
Grade	RUO - IVD	RUO - IVD	RUO - IVD	RUO - IVD