

Instruments for Western Blotting

BlueVertical™ PRiME™ & TankBlotter


- Dual mini tank systems
- Accommodates 1 – 2 gels in cassettes with outer dimensions of 10 cm x 10 cm x 0.7 cm
- Leak-free inner core unit with a unique clamp system
- Unique, easy-to-handle blotting insert for tank blotting
- Gel Casting Stand for pouring two gels
- Quality designed and made in Germany



Product	Size	Cat. no.
BlueVertical™ PRiME™	1 unit	BV-104
BlueVertical™ PRiME™ TankBlot	1 unit	BV-104-TB
BlueVertical™ PRiME™ Casting Stand	1 unit	BV-104-CS
BlueVertical™ PRiME™ Blot Module	1 unit	BV-104-B

BlueBlot Semi-Dry Blotter

- Best blot results for small and large proteins
- Platinum-covered steel net as anode, spring-mounted
- Stainless steel plate as cathode
- Easy to clean
- Versatile models:
 - 11 cm x 11 cm for mini gels
 - 17 cm x 17 cm for mini and midi gels
 - 26 cm x 24 cm for 2D HPE™ BlotGels and other large format gels
- Deployable for thicker gels and blotting stacks



Product	Blot Area	Size	Cat. no.
BlueBlot Semi-Dry Blotter SD 11	11 cm x 11 cm	1 unit	BB-SD11
BlueBlot Semi-Dry Blotter SD 17	17 cm x 17 cm	1 unit	BB-SD17
BlueBlot Semi-Dry Blotter SD 26	24 cm x 26 cm	1 unit	BB-SD26
Electrode Set for BB-SD11	11 cm x 11 cm	1 unit	BB-E11
Electrode Set for BB-SD17	17 cm x 17 cm	1 unit	BB-E17

BluePower™ Power Supply for Blotting

- 300 V, 2000 mA, 300 W
- Semi dry blotting applications including large format gels up to 24 cm x 26 cm
- Fully programmable (9 programmes with 9 steps each)
- 4x 2 outlets



Product	Size	Cat. no.
BluePower™ 300 BLOT Power Supply	1 unit	BP-300-BLO